

Details for Users

Ordering Information

For ordering CBO128 please give **Article Number**, **Cap Cut**, and **Size**

- Article Number: **CBOL128**
- Cap Cut: **Caucasian** or **Asian**
- Size (given in cm head circumference):
 - Adult caps: **54, 56, 58, 60, 62, 64** (average male: 58, average female: 54)
 - Children caps: **52** (7 years), **54** (11 years)

The catalogue-number comprises the cap as described, serial number, and this document; all packed in a labelled cardboard box. For further information about accessories or consumables, please visit our website or contact our local distributor.

Cap

Sublunion Cap with integrated chin belt, high-precision fabric, black, occipital velcro patch.

Options: Caucasian or Asian, Size (for additional caps)

Cap is equipped with A07 electrode holders, with CGX labels according to the layout.

Cleaning and Maintenance

Please refer to the Mobile-128 manual

Coordinates of Electrode Positions

Please contact info@easycap.de to receive Theta/Phi- or XYZ-coordinates.