



## Details for Users

### Ordering Information

For ordering CBOL72 please give **Article Number**, **Cap Cut**, and **Size**:

- Article Number: **CBOL72**
- Cap Cut: **Caucasian** or **Asian**
- Size (given in cm head circumference):
  - Adult caps: **54, 56, 58, 60, 62, 64** (average male: 58, average female: 54)
  - Children caps: **50** (5 years), **52** (7 years), **54** (11 years)
  - Infant caps: **42, 44** (7 month), **46, 48** (3 years)

The catalogue-number comprises the cap as described, serial number, and this document; all packed in a labelled cardboard box. For further information about accessories or consumables, please visit our website or contact our local distributor.

### Cap

SubInion Cap with integrated chin belt, high-precision fabric, black, occipital velcro patch.

Options: *Caucasian* or *Asian*, *Size*

Cap is equipped with A07 electrode holders, with CGX labels according to the layout.

### Cleaning and Maintenance

Please refer to the Mobile-72 manual.

### Coordinates of Electrode Positions

Please contact [info@easycap.de](mailto:info@easycap.de) to receive Theta/Phi- or XYZ-coordinates.