

actiCAP Xpress Twist

dry, fast & comfortable



- ✓ Active dry electrode system for fast and zero-gel recordings
- ✓ Individual electrode adjustment for maximum comfort
- ✓ Gold-coated electrode tips in different lengths and shapes for high signal quality
- ✓ Cap sizes range from infants to adults; available channel counts: 8, 16, 32



Unique twist mechanism

- ✓ Easy to use twisting mechanism built into the electrode holders allows you to individually adjust the sensor's height
- ✓ Optimal connection to the participant's scalp
- ✓ Achieve high-quality data and maintain participant comfort

High quality dry electrodes

- ✓ Active dry electrodes based on more than 15 years of experience with active electrode technology
- ✓ Gold plated exchangeable dry sensors for lowest noise levels)

Soft and flexible caps

- ✓ 8, 16 or 32 channel standard cap layouts
- ✓ Caps are made of neoprene and can be easily customized by the user with a total of 78 pre-defined electrode positions
- ✓ Additional cutouts reduce sweat artifacts and allow access to sensors even when the cap is already mounted

