



BRAIN PRODUCTS

Solutions for neurophysiological research

BrainAmp dualBUA Supplement

Brain Products GmbH
Zeppelinstrasse 7
82205 Gilching
Germany

T +49 (0) 8105 733 84 0
F +49 (0) 8105 733 84 505
www.brainproducts.com

Worldwide Caution:

Our products are scientific equipment
for INVESTIGATIONAL USE ONLY!

Medical use e.g. for diagnosis, treatment
of disease or other such purposes is
strictly forbidden.





Use the dualBUA adapter cable to connect the two USB2 Adapters* to one another as shown in the figure. Make sure you also tighten the retaining screws on the connectors. The dualBUA is intended to be used only for combinations using 2 x BUA128 or 1 x BUA128 and 1 x BUA64.

**Compatible with BrainVision Recorder Version 1.20.0401 or higher.*

