


## actiCAP 128Ch Standard-2

- Green holders: electrode positions Ch1 – Ch32
- Yellow holders: Label 1-32, hard-wired Ch33 – Ch64
- Red holders: Label 1-32, hard-wired Ch65 – Ch96
- White holders: Label 1-32, hard-wired Ch97–Ch128
- Blue holder: Ref
- Black holder: Gnd

## Components:

- Softcap, white
- Holders with flat side inwards
- 3 additional holders (white, empty label) for use with double-sided adhesive rings to place electrodes on bare skin
- 1 chin belt
- this layout / pinout



Contents of these documents are the intellectual property of Brain Products GmbH, and are subject to change without specific notification. Brain Products GmbH does not grant warranty or assume liability for the correctness of individual statements herein. Nor does Brain Products GmbH enter into any obligation with regard to these documents.